

Protocol for Preparing 1M IPTG Stock Solution (0.5 mL Aliquots)

Materials:

- IPTG (Isopropyl β -D-1-thiogalactopyranoside) powder
 - Deionized water (DI water)
 - Sterile 1.5 mL microcentrifuge tubes
 - 0.22 μ m sterile syringe filter or 0.22 μ m sterile filter unit
 - Pipettes and sterile tips
 - Scale for weighing
-

Steps:

1. **Weigh the IPTG:**
 - Measure **2.38 g** of IPTG powder (molecular weight: 238.3 g/mol) using a clean, calibrated scale. This amount will make 10 mL of a 1M solution.
 2. **Dissolve in deionized water:**
 - Transfer the IPTG powder into a sterile container.
 - Add **10 mL** of deionized water.
 - Gently mix by pipetting up and down or using a vortex until the IPTG is fully dissolved.
 3. **Sterilize the solution:**
 - Use a **0.22 μ m sterile filter** to sterilize the IPTG solution. Draw the solution into a syringe, attach the filter, and expel the solution into a sterile container.
 4. **Aliquot the solution:**
 - Pipette **0.5 mL** of the sterilized 1M IPTG solution into each sterile 1.5 mL microcentrifuge tube.
 5. **Label the aliquots:**
 - Clearly label each tube with "1M IPTG," the date, and your initials.
 6. **Storage:**
 - Store the aliquots at **-20°C** for long-term use. Aliquots can be thawed as needed.
-

Note:

- IPTG is stable for several months at -20°C. Avoid repeated freeze-thaw cycles by only thawing the aliquot needed for use.